

Limb bud micromass culture

Karen M. Lyons Lab – UCLA

Previously described in H.-S. Seo. R. Serra, Dev Biol, 2007, 310(2):304-316

Modified on 04/05/10

Reagents/solutions:

- 1) Collagenase D, Roche Applied Sci, Cat# 11088858001
- 2) DPBS (1x), Mediatech, Inc., Cat# 21-031-CV
- 3) F12 (Ham, 1x), Invitrogen, Cat# 11765062
- 4) DMEM media, Mediatech, Inc., Cat# 10-017-CM
- 5) β -Glycerol phosphate (Sigma, #G9891)
 - a. To make 10 ml of 1M β -Glycerol phosphate
 - i. Add 2.1604 g of β -Glycerol phosphate to 10 ml DPBS
 - ii. Filter-sterilize
 - iii. Store at 4 °C
- 6) Ascorbic acid (Sigma, #A4544)
 - a. To make 50 mg/ml of ascorbic acid
 - i. Add 500 mg of ascorbic acid to 10 ml of DMEM media
 - ii. Filter-sterilize
 - iii. Store 1 ml aliquots (covered in foil) at -20 °C

Procedure:

- 1) Isolate limb buds from E11.5 day mouse embryos
- 2) Dissociate mesenchymal cells into single cell suspension by digesting in 300-500 μ l of 1 mg/ml collagenase D (diluted from 100 mg/ml stock in sterile DPBS) for 1-2 hour in 37 °C incubator
 - a. Vortex every 30 min to facilitate cell dissociation
- 3) Filter through 70 μ m filter
- 4) Pellet cells by centrifugation at 1000 rpm for 5 min
- 5) Aspirate collagenase solution
- 6) Resuspend in 50 μ l F12:DMEM (1:1) media supplemented with 10 % FBS and antibiotics
- 7) Reconstitute at density of 1×10^7 cells/ml
- 8) Drop 10 μ l of cell suspension into each well of a 24-well plate
- 9) Allow cells to attach to dish by incubating for 1 hr in 37 °C incubator
- 10) After incubation, flood cells with 500 μ l of chondrogenic media (F12:DMEM (1:1) media supplemented with 10 % FBS and antibiotics, 50 μ g/ml ascorbic acid, 10 mM β -Glycerol phosphate) **supplemented with appropriate growth factors, inhibitors, etc.**
- 11) Incubate cultures in 37 °C incubator for up to three days
 - a. Can process cells for Alcian blue staining, immunofluorescence staining, western blotting, or real-time PCR